

Emerging Contaminants

Comprehensive Targeted Workflows

Chris Snelling
Agilent Technologies




PFAS

Per/Poly Fluoroalkyl Substances


Terminology

Perfluoroalkyl substance



ALL H atoms linked to C in alkyl chain are substituted with F

Polyfluoroalkyl substance




SOME (*but not all*) H atoms linked to C in alkyl chain are substituted with F

PFAS Classifications and Terminology

>4000 PFAS compounds in commerce

Common Acronyms

PFCA	Perfluoroalkylcarboxylic acid
PFOA	Perfluorooctanecarboxylic acid
PFAS	Perfluoroalkylsulfonate
PFOS	Perfluorooctanesulfonate
PFASi	Perfluoroalkylsulfinate
FOSA	Perfluorooctanesulfonamide
FOSAA	Perfluorooctanesulfonamidoacetic acid
FOSE	Perfluorooctanesulfonamidoethanol
FTOH	Fluorinated telomer alcohol (-OH functional group)
FTA	Fluorinated telomer acid
FTUA	Fluorinated telomer unsaturated acid
FTS	Fluorinated telomer sulfonate
PFAPA	Perfluoroalkylphosphonic acid
PFPI	Perfluoroalkylphosphinate
PAP	Mono-substituted polyfluoroalkylphosphate ester
diPAP	Di-substituted polyfluoroalkylphosphate ester
PFAI	Perfluoroalkyl iodide
SFA	Semifluorinated alkane
FTI	Fluorinated telomer iodide
FTO	Fluorinated telomer olefin
FTAC	Fluorinated telomer acrylate



Wang, Z et al. (2017). *Environ. Sci. Technol.* 51, 2508-2518.

Current Standard & Consensus Methods Available

Method	Matrix Tested	No. of analytes	Sample preparation procedure	Quantification Technique	Sample (mL)/Injection (uL) Volume	Year
EPA 537	Drinking Water	14	Solid phase extraction	Internal standard correction	250/10	2008
EPA 537.1	Drinking Water	18	Solid phase Extraction	Internal standard correction	250/10	2018
EPA 533	Drinking Water	25	Solid phase extraction	Isotope dilution	250/10	2019
EPA 8327	Surface water, Ground water, Wastewater influent and effluent	24	Dilute & shoot	External calibration (isotope dilution also allowed)	5/30	2019
EPA 1633 (3 rd Draft)	Wastewater, Soil, Biota, Sediment, Groundwater	40	Solid Phase Extraction	Isotope Dilution	500/2	2022
ASTM 8421	Surface water, Ground water, Wastewater influent and effluent	44	Dilute & shoot	External calibration (isotope dilution also allowed)	5/30	2022
ASTM 7968	Soil and solids	21	Organic extraction with MeOH	External calibration	5g/30	2015
ISO/DIS 21675:2019	Drinking Water, Sea water, Fresh water, wastewater (<0.2% solids)	30	Solid phase extraction	Internal standard correction	500/5	2019

National Primary Drinking Water Regulations

Selected Per- and poly-fluoroalkyl substances (PFAS)

- [PFAS rule references](#)

Contaminant	MCLG ¹ (mg/L) ²	MCL ¹ (mg/L) ²	HBWC ⁹ (mg/L) ² for Hazard Index Calculation
Hazard Index PFAS (HFPO-DA, PFBS, PFHxS, and PFNA)	1(unitless)	1(unitless)	Not applicable
HFPO-DA (commonly known as GenX Chemicals)	0.00001	0.00001	0.00001
PFBS	No individual MCLG	No individual MCL	0.002
PFHxS	0.00001	0.00001	0.00001
PFNA	0.00001	0.00001	0.00001
PFOA	zero	0.0000040	Not applicable
PFOS	zero	0.0000040	Not applicable

4 ng/L! (4 ppt)



TNI PT for Accreditation
Fields of Proficiency Testing with PTRLs
Drinking Water
Effective: January 1, 2025

Blue = New Analyte

Magenta = Changes

Matrix	EPA Analyte Code	TNI Analyte Code	CAS Number	Analyte ²	Conc Range	Acceptance Criteria ^{3,4,5,6}				TNI PTRL ⁷
						a	b	c	d	
				PAH	µg/L					µg/L
Drinking Water	0122	5580	50-32-8	Benzo(a)pyrene ¹	0.2 to 2.5	0.8471	-0.0040	0.1854	0.0547	0.02
				Dioxin	pg/L					pg/L
Drinking Water	0252	9618	1746-01-6	2,3,7,8-Tetrachlorodibenzo- p-dioxin (2,3,7,8-TCDD) ¹	20 to 100	0.8642	1.4865	0.1392	1.1445	11
				PFAS	ng/L					ng/L
Drinking Water	2813	9490	763051-92-9	11-Chloroeicosafuoro-3-oxaundecane-1-sulfonic acid (11Cl-PF3OUdS)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2822	6948	39108-34-4	1H,1H,2H,2H-Perfluorodecane sulfonic acid (8:2FTS)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2821	6946	757124-72-4	1H,1H,2H,2H-Perfluorohexane sulfonic acid (4:2FTS)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2820	6947	27619-97-2	1H,1H,2H,2H-Perfluorooctane sulfonic acid (6:2FTS)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2815	6951	919005-14-4	4,8-Dioxo-3H-perfluorononanoic acid (ADONA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2814	6952	756426-58-1	9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid (9Cl-PF3ONS)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2816	9460	13252-13-6	Hexafluoropropylene oxide dimer acid (HFPO-DA) (GenX)	10 to 200			±30% fixed acceptance limit		7
Drinking Water	2817	4846	2991-50-6	N-ethyl perfluorooctanesulfonamidoacetic acid (NEtFOSAA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2818	4847	2355-31-9	N-methyl perfluorooctanesulfonamidoacetic acid (NMeFOSAA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2827	6956	151772-58-6	Nonafluoro-3,6-dioxaheptanoic acid (NFDHA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2826	6957	113507-82-7	Perfluoro(2-ethoxyethane)sulfonic acid (PFEESA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2823	6965	377-73-1	Perfluoro-3-methoxypropanoic acid (PFMPA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2825	6966	863090-89-5	Perfluoro-4-methoxybutanoic acid (PFMBA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2801	6918	375-73-5	Perfluorobutanesulfonic acid (PFBS)	10 to 200			±30% fixed acceptance limit		7
Drinking Water	2819	6915	375-22-4	Perfluorobutanoic acid (PFBA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2807	6905	335-76-2	Perfluorodecanoic acid (PFDA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2808	6903	307-55-1	Perfluorododecanoic acid (PFDoA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2829	9470	375-92-8	Perfluoroheptanesulfonic acid (PFHpS)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2802	6908	375-85-9	Perfluoroheptanoic acid (PFHpA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2803	6927	355-46-4	Perfluorohexanesulfonic acid (PFHxS)	10 to 200			±30% fixed acceptance limit		7
Drinking Water	2809	6913	307-24-4	Perfluorohexanoic acid (PFHxA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2804	6906	375-95-1	Perfluorononanoic acid (PFNA)	10 to 200			±30% fixed acceptance limit		7
Drinking Water	2805	6931	1763-23-1	Perfluorooctanesulfonic acid (PFOS)	10 to 200			±30% fixed acceptance limit		7
Drinking Water	2806	6912	335-67-1	Perfluorooctanoic acid (PFOA)	10 to 200			±30% fixed acceptance limit		7
Drinking Water	2828	6934	2706-91-4	Perfluoropentanesulfonic acid (PFPeS)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2824	6914	2706-90-3	Perfluoropentanoic acid (PFPeA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2810	6902	376-06-7	Perfluorotetradecanoic acid (PFTTA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2811	9563	72629-94-8	Perfluorotridecanoic acid (PFTrDA)	10 to 200			±40% fixed acceptance limit		6
Drinking Water	2812	6904	2058-94-8	Perfluoroundecanoic acid (PFUnA)	10 to 200			±40% fixed acceptance limit		6



Rulemaking at DEQ

DEQ's mission is to be a leader in restoring, maintaining and enhancing the quality of Oregon's air, land and water.

› [Rulemaking at DEQ](#) › [Proposed Rules](#) › [PFAS 2025](#)

PFAS 2025

How to Find Rules and Rulemaking Documents

Search DEQ Rulemaking Documents

Search Official Version Of The Rules

Greenhouse Gas Emissions Program 2021

Proposed rule


Rulemaking contact: [Sarah Van Glubt](#), 503-709-8253, PFAS.2025@deq.oregon.gov

This rulemaking proposes to include two per- and polyfluoroalkyl substances (PFAS), perfluorooctanoic acid (PFOA) and perfluorooctanesulfonic acid (PFOS), including their salts and structural isomers, in the definition of hazardous substances in Oregon Administrative Rule (OAR) 340-122-0115 (30) which would give the Department authority to require investigation and removal and remedial actions of PFOA and PFOS releases and align with the US Environmental Protection Agency's approach. PFOA and PFOS are among the most commonly detected PFAS, are known or suspected carcinogens, and have been shown to exhibit toxicity effects to humans and wildlife even at low levels of exposure.

PFAS Sample Prep & Analysis: Basic Workflow Overview

Address background and contamination issues with a robust workflow

PFAS LC/MS Workflow



SLIMS - Automation of EPA Workflows for PFAS

Sample Management

Inventory Management

Workflow Execution

- Sample Prep
- Masshunter Integration
- Results Evaluation

Dashboard & Reporting



Workflow Graph

Protocol: Soil Sample Prep Version: 1

Seq...	Name	Type	Plugin	Renderer	Descrip...	Tests	Conten...	Sub pr...	Sign off...
1	Link Sample Prep Kit	Link content			Link and consume				
2	Adjust pH	Result			pH should be 6.0 - 7.0	pH			
3	Extraction Setup	Mix			Clean silanized		Output		
4	Condition SPE	Mix			15 mL - 1% methanolic		Output		

Sort: Seq no (Ascending) Total rows: 5


Protocol step info

Name	Dr.	Usage t...	Default ...	Defe
InfinityLab Ultrapure LC/MS Water		<input checked="" type="checkbox"/>	500	ml
InfinityLab Ultrapure LCMS Acetonitrile		<input checked="" type="checkbox"/>	250	ml
LC-MS grade MeOH		<input checked="" type="checkbox"/>	25	ml


PFAS Analysis Workflow Consumables, Supplies and Hardware



Wastewater
(500mL)



Centrifuge tubes
50 mL (5610-2039)



Carbon S
(5610-2093)



Empty SPE tubes,
60 mL (12131012)



SPE Adapters
(12131001)



Glass Wool,
Silanized
(8500-1572)



Bond Elut PFAS
WAX 150 mg, 6 mL
(5610-2152)



Vac Elut SPS 24
(12234003)




PFC-free kit (5004-0006)
PFC delay column (5062-8100)
Analytical column Eclipse Plus
C18, 2.1 x 100 mm, 1.8 μm
(959758-902)

PFAS MRM Database (G1736AA)




1290 Infinity II LC
system

6495 triple
quadrupole LC/MS



Polypropylene AS Vials
(5191-8151, 5191-8150)



Syringes
(9301-6476)



Nylon Filters
(5190-5092)



Centrifuge tubes
15 mL (5610-2039)

Solid Phase Extraction – Bond Elut PFAS WAX

Specifically designed, developed and manufactured for PFAS applications



- Cleanliness
- Sorbent and cartridge formats compatible with all existing regulated methods
 - EPA method 533 for drinking water
 - EPA method 1633 (draft) for aqueous, solids, biosolids, and tissue samples
 - ISO 21675:2019 for drinking water, sea water, fresh water, and wastewater
- Performance equivalent to other commercial cartridges
- Fits into Agilent's existing PFAS workflows



New way to overcome strong solvent effects at high injection volumes

Improve peak shape and sensitivity of early eluting PFAS

- New 1260 Infinity II Hybrid Multisampler (G7167C)
- Feed injection valve
- Allows dilution during injection
- Classic injection mode simulation by fast feed speed




The 6495 Triple Quadrupole LC/MS with iFunnel technology


Reliability & Robustness for demanding applications.

Born out of our drive for innovation and improvement, Agilent's 6495 Triple Quadrupole LC/MS is an innovative, powerful, and reliable mass spectrometer for customer applications demanding the lowest limits of detection, sampling speed, and signal reproducibility.

Experience the highest level of confidence with the 6495 Triple Quadrupole LC/MS.




Comprehensive Workflow for PFAS Analysis in Wastewater




Poster Thursday 075

Method Detection Limit and Relative Standard Error




6495 LC/TQ Robustness

>300 injections run over 80+ hours; RSD <6% for all 40 PFAS in EPA 1633



EPA 533




	3 rd Generation iFunnel LOQ (ng/L)	4 th Generation iFunnel LOQ (ng/L)
PFBS	0.01	0.004
HFPO-DA	0.01	0.01
PFOA	0.02	0.01
PFHxS	0.02	0.01
PFNA	0.03	0.03
PFOS	0.2	0.03


2.5- 7x Improvement

PFAS Analysis (ASTM 8421)

- Dilute and Shoot method
- Reporting range 10 – 400 ng/L
- Injection volume: **10 uL**



Enhanced Sensitivity with 6495!




PFAS Workflows Conclusions


- A comprehensive workflow including sample management, preparation, consumables, data acquisition/analysis, and reporting developed for the PFAS analysis.
- Reliable sample preparation with excellent recovery.
- 6495 LC/ TQ showed excellent reproducibility and robustness for targeted analysis
- Instrument sensitivity allowed the reduced injection volume for dilute and shoot method
- **Now a Fully Mature Workflow**

6-PPD Quinone

p-Phenylenediamine Quinone




Oxidative Formation of 6-PPD Quinone from 6-PPD



Species	LC ₅₀ (µg/L)	Test duration (h)	Toxicity Key
Coho salmon (<i>Oncorhynchus kisutch</i>)	0.04, ²⁴ 0.08, ²⁵ 0.095 ²	24	Higher
White-spotted char (<i>Salvelinus leucomaenis pluvius</i>)	0.51 ²⁶	24	
Brook trout (<i>Salvelinus fontinalis</i>)	0.59 ³	24	
Rainbow trout/steelhead (<i>Oncorhynchus mykiss</i>)	0.64, ²⁹ 1.0, ³ 2.26 ⁵	96	
Chinook salmon (<i>Oncorhynchus tshawytscha</i>)	67.3 ²⁴ , 82.1 ²⁵	24	
Sockeye salmon (<i>Oncorhynchus nerka</i>)	Not acutely toxic at 50 ²⁵	24	Lower
Atlantic salmon (<i>Salmo salar</i>)	Not acutely toxic at 12.2 ²⁸	48	
Brown trout (<i>Salmo trutta</i>)	Not acutely toxic at 12.2 ²⁸	48	
Arctic char (<i>Salvelinus alpinus</i>)	Not acutely toxic at 12.7 ³	24	
Southern Dolly Varden (<i>Salvelinus curilus</i>)	Not acutely toxic at 3.8 ²⁶	48	
Cherry salmon (<i>Oncorhynchus masou masou</i>)	Not acutely toxic at 3.5 ²⁶	48	

DRAFT Method 1634

Determination of 6PPD-Quinone in Aqueous Matrices Using Liquid Chromatography with Tandem Mass Spectrometry (LC/MS/MS)



6-PPD Quinone in Stream Water Method via LC-MS/MS

Compound Name	Precursor Ion (m/z)	Product Ion (m/z)	Ret. Time (min)	Fragmentor (V)	Collision Energy (V)	Cell Accelerator Voltage	Polarity
6PPD-Quinone	299.2	241.1	3.52	105	32	4	Positive
6PPD-Quinone	299.2	215.1	3.52	105	16	5	Positive
6PPD-Quinone	299.2	187.1	3.52	105	32	5	Positive
D ₅ -6PPD-Quinone	304.2	246.1	3.49	110	36	4	Positive
D ₅ -6PPD-Quinone	304.2	220.1	3.49	110	20	4	Positive
D ₅ -6PPD-Quinone	304.2	192.1	3.49	110	36	5	Positive





Figure 1. 6PPD-quinone chromatography at 0.2 ng/mL.

6PPD-quinone was quantitated using a calibration curve from 0.01 to 50 ng/mL using quadratic fit, 1/x weighting, and including the origin. The R² value was greater than 0.999.



Method recovery (accuracy) in HPLC water was 113.5% and in stream water was 112.6%. Precision, expressed as %RSD, was 3% in HPLC water and 1% in stream water for n = 5 replicates at 5 ng/mL. The LCMRL of 6PPD-quinone was 0.023 ng/mL. There was no contribution to the 6PPD-quinone signal from the LC system.

6-PPD Quinone In Salmon Tissue



EMR sorbent technology effectively traps lipids through two mechanisms:

- ✓ **Size exclusion** – Unbranched hydrocarbon chains (lipids) enter the sorbent; bulky analytes do not
- ✓ **Sorbent chemistry** – Lipid chains that enter the sorbent are trapped by hydrophobic interactions



LC-MS/MS

New Contaminants, New Tools, New Solutions





Agilent

Trusted Answers